
Accessing the Taltio AIC Network

1. On the microscope


a. Files can be saved directly in Taltio AIC instrument drive or moved there after the imaging session.



b. If the network location is not there, click on the connect icon on the desktop.




c. Make sure to create a folder for your research group to keep files organized for everyone.




Do not leave data in the Taltio AIC Network!
FILES OLDER THAN 30 DAYS WILL BE AUTOMATICALLY DELETED!!

2. Accessing data on AIC processing workstation

a. You can access all data in Taltio AIC Network from the processing workstations Alpha and Beta for image processing and analysis. On the workstations, You'll find an AIC network drive




b. Within the AIC drive, each instrument has an assigned folder, with all data synced




3. Accessing files from your university computer

a. Access to Taltio is available to University of Turku researchers. For Åbo Akademi University researchers, please reach out to IT Helpdesk to possibly request a Taltio visitor account.


b. Access OpenIris and select the Resources tab



c. Search for “Taltio AIC Network | Access request”




d. Enter your email, submit your request, and wait for request confirmation




e. Once enabled as a user, in windows File explorer, click on This PC




f. Click on the three dots menu




g. Select “Map network drive”




h. As folder enter: \\utu.fi\taltio\Instruments\bioscience\AIC, then click Finish



i. You should now see an AIC drive in your PC



j. In the AIC, each microscope has its own folder



k. Transfer data from the Taltio AIC Network to your computer or Lab Taltio folder or other storage location.
For security reasons, users can only read files from their personal computer but not modify or delete them.